EN

DAMALON

Polyamide 6-6.6 (Nylon)
Monofilament, Multifilament, Synthetic
Non-absorbable sterile surgical suture
Black, Blue, Undyed

DESCRIPTION
DAMALON suture is synthetic non-absorbable sterile surgical suture. It is composed of Polyamide polymers. It is Monofilament Multifilament dyed sutures. **DAMALON** suture has high in vivo tensile strength, does not support bacterial growth and has the further advantages of being remarkably smooth & easy to handle. **DAMALON** meets all requirements for non-absorbable synthetic sutures defined by the European Pharmacopoeia, EC Medical Device Regulation 2017/745/AT, and US Pharmacopoeia. Safety and clinical performance are detailed in GMD.SSCP.07.

INDICATIONS
DAMALON is indicated for general soft tissue approximation and/or ligation (including gastrointestinal, obstetric and gynecological and urological surgery), including use in ophthalmic surgery, and for the closure and closure of suturable wounds by initial wound suturing and stapling under sterile conditions.

CONTRAINDICATIONS
DAMALON suture should not be used where permanent retention of tensile strength is required as in fixation of intraocular lenses or synthetic vascular grafts, due to gradual loss of tensile strength which may occur over prolonged periods in vivo.

WARNINGS
Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing **DAMALON** suture for wound closure, as risk of application and the suture material used. As with any foreign body, prolonged contact of any suture with soft tissues, such as those found in the urinary or biliary tracts may result in calculus formation. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. Discard opened packages and unused sutures.
Do not re-use.
Don not re-sterilize.

PRECAUTIONS
In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps of needle holders. **DAMALON** sutures, which are treated to enhance handling characteristics, require the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of addition throws may be appropriate when knotting monofilament polyamide sutures.

To avoid damaging needle points and swage areas, grasp the needle in the area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.
Discard used needles in "sharps" container. Avoid prolonged exposure to elevated temperatures.
Don't use after expiry date.

ADVERSE REACTIONS
Adverse effects associated with the use of this device include wound dehiscence, minimal acute inflammatory tissue reaction, calculus formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile accours, gradual loss of tensile strength over time, enhanced bacterial infectivity pain, and transitory local irritation and edema at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needlesticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

STERILITY
DAMALON sutures are sterilized by ethylene oxide. Sterility is preserved only when opened under sterile conditions. Do not re-sterilize.
Do not use if package is opened or damaged.
Discard opened unused sutures.

STORAGE
Keep away from moisture and direct heat. Recommended storage condition is between 5-25°C.
Don't use after expiry date.

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