EN DAMAESTER

Polyester
Braided, Coated
Non-absorbable sterile surgical suture

DESCRIPTION
DAMAESTER suture is non-absorbable sterile surgical suture, composed of fine filaments of polyethylene terephthalate (polyester). It is braided, silicone coated, dyed or undyed sutures. **DAMAESTER** suture has excellent strength properties, and it is Coated with Silicone to provide for smoother and better handling. **DAMAESTER** meets all requirements for non-absorbable synthetic sutures defined by the European Pharmacopoeia, EC Medical Device Regulation 2017/745/AT, and US Pharmacopoeia. Safety and clinical performance are detailed in GMD.SSCP.08

INDICATIONS
DAMAESTER sutures are indicated for use in soft tissue approximation and/or ligation, including ophthalmic procedures. **DAMAESTER** sutures are also indicated for use in cardiovascular surgery.

CONTRAINDICATIONS
DAMAESTER suture, should not be used where permanent retention of tensile strength is required as in fixation of intraocular lenses or synthetic vascular grafts, due to gradual loss of tensile strength which may occur over prolonged periods in vivo.

WARNINGS
Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing **DAMAESTER** suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. Discard opened packages and unused sutures. Do not re-use. Do not re-sterilize.

PRECAUTIONS
In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. **DAMAESTER** sutures, which are treated to enhance handling characteristics, require the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be appropriate when knotting coated polyamide sutures. To avoid damaging needle points and swage areas, grasp the needle in the area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point.

Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" container. Avoid prolonged exposure to elevated temperatures. Don't use after expiry date.

ADVERSE REACTIONS
Adverse effects associated with the use of this device include wound dehiscence, minimal acute inflammatory tissue reaction, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, transitory local irritation at the wound site and all foreign bodies it may potentiate an existing infection. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needlesticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

STERILITY
DAMAESTER sutures are sterilized by ethylene oxide. Sterility is preserved only when opened under sterile conditions. Do not re-sterilize. Do not use if package is opened or damaged. Discard opened unused sutures.

STORAGE
Keep away from moisture and direct heat. Recommended storage condition is between 5-25°C. Don't use after expiry date.

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