EN

DAMADIOXANONE

Poly (P-dioxanone)
Monofilament Synthetic
Absorbable sterile surgical suture

DESCRIPTION
DAMADIOXANONE suture is synthetic absorbable sterile surgical suture. It is composed of polymers made from 100% homo polymer of polydioxanone. It is Monofilament sutures dyed (violet) or undyed. **DAMADIOXANONE** suture is non-antigenic, non-pyrogenic and elicits only a mild tissue reaction during the absorption process, and is completely absorbed in 180-240 days. The in vitro retention of strength is 65% in four weeks. **DAMADIOXANONE** meets all requirements for absorbable synthetic sutures defined by the European Pharmacopoeia, EC Medical Device Regulation 2017/745/AT, and US Pharmacopoeia. Safety and clinical performance are detailed in GMD.SSCP04

INDICATIONS
Polydioxanone monofilament synthetic absorbable sutures are indicated for use in all types of soft tissue approximation (including gastrointestinal, obstetric, gynecologic, and urologic surgery), including tissue where growth is anticipated, and in ophthalmic surgery. Polydioxanone sutures are not indicated for adult cardiovascular tissue, microsurgery, or neural tissue. These sutures are particularly useful when the combination of an absorbable suture and extended wound support (up to six weeks) is desired.

CONTRAINDICATIONS
DAMADIOXANONE suture, being absorbable, should not be used where prolonged approximation of tissue under stress is required (beyond six weeks), or in conjunction with prosthetic devices, i.e., heart valves or synthetic grafts.

WARNINGS
The safety and effectiveness of Polydioxanone have not been established in neural tissue, adult cardiovascular tissue or for use in microsurgery. Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing **DAMADIOXANONE** suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Users should consider the in vivo performance when selecting a suture. The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, callus formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site.

Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needlesticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

STERILITY
DAMADIOXANONE sutures are sterilized by ethylene oxide. Sterility is preserved only when opened under sterile conditions. Do not re-sterilize. Do not use if package is opened or damaged. Discard opened unused sutures.

As with any foreign body, prolonged contact of any suture with soft tissues, such as those found in the urinary or biliary tracts may result in calculus formation. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. Discard opened packages and unused sutures. Do not re-use. Do not re-sterilize.

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STORAGE
Keep away from moisture and direct heat. Recommended storage condition is between 5-25°C. Don't use after expiry date.

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