EN
DAMASILK
Non-absorbable Silk Braided, Coated Sterile Surgical Suture

DESCRIPTION
DAMASILK suture is a natural non-absorbable sterile surgical suture. It is composed of approximately 70% fibrin protein fibers and 30% extraneous material (gum). A specially developed process removes extraneous material amounting to about 30% of the original weight of the raw silk. It is a braided, silicone-coated, black-dyed suture. **DAMASILK** suture elicits only a minimal initial tissue reaction followed by gradual encapsulation by fibrous connective tissue, and loses most of its tensile strength within one year in vivo. **DAMASILK** suture is not absorbed nor is it subjected to degradation or weakening by the action of tissue enzymes. It is coated with silicone to allow easy passage through tissue. **DAMASILK** meets all requirements for non-absorbable synthetic sutures defined by the European Pharmacopoeia, EC Medical Device Regulation 2017/745/AT, and US Pharmacopoeia. Safety and clinical performance are detailed in GMD.SSCP.06

PRECAUTIONS
In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

ADVERSE REACTIONS
Adverse reactions associated with the use of this device include wound dehiscence; gradual loss of tensile strength over time; allergic response in patients who are known to be sensitive to silk; calculus formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs; transitory local irritation at the wound site; and like all foreign bodies it may potentiate an existing infection. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needlesticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

INDICATIONS
DAMASILK is indicated for general soft tissue approximation and/or ligation (including gastrointestinal, obstetric and gynecological and urological surgery), closure and dressing of suturable wounds under sterile conditions by initial suturing of the wound.

CONTRAINDICATIONS
DAMASILK is contraindicated in patients with known allergies or sensitivities to silk. It should also not be used where permanent retention of tensile strength is required, such as in fixation of intraocular lenses or with synthetic vascular grafts, due to gradual loss of tensile strength that may occur over prolonged periods in vivo. Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing **DAMASILK** suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. As with any foreign body, prolonged contact of any suture with salt solutions such as those found in the urinary or biliary tracts may result in calculus formation. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. Discard opened packages and unused sutures. Do not re-use. Do not re-sterilize.

WARNINGS
DAMASILK sutures, which are treated to enhance handling characteristics, require the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and the experience of the surgeon.

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