DAMACRYL VET

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DAMACRYL VET
GE000012 REV-10
Date of issue: 12/2018
Date of revision: 07/2025

DAMACRYL VET

Polyglycolic Acid
Coated & Braided (Violet) and Undyed
Synthetic absorbable sterile surgical suture

DESCRIPTION

DAMACRYL VET suture is synthetic absorbable sterile surgical suture. It is composed of homo polymer prepared and synthesized from 100% Polyglycolic Acid. Thess sutures are braided, dyed (violet) or undyed and coated with polycaprolactone and calcium stearate.

DAMACRYL VET is non-antigenic and non-pyrogenic.

During absorption, only minimal tissue reaction occurs and it is fully absorbed by the body within 60-90 days. In vitro, it retains 65% of its strength during the first 2 weeks.

DAMACRYL VET meets all requirements for absorbable synthetic sutures defined by the European Pharmacopoeia, EC. Medical Device Regulation 2017/745, and US Pharmacopoeia. Safety and clinical performance are detailed in GMD.SSCP.01.

INDICATIONS

DAMACRYL VET and PGA are used for approximation or ligation of general soft tissues (including gastrointestinal, obstetric, gynecologic, and urologic surgeries), closure of wounds under sterile conditions, and stapled or sutured wounds.

DAMACRYL VET is not indicated for cardiovascular or neurological tissues.

CONTRAINDICATIONS

DAMACRYL VET suture, being absorbable, should not be used where extended approximation of tissue is required.

WARNINGS

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing **DAMACRYL VET** suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Users should consider the in vivo performance when selecting a suture. The use of this suture may be inappropriate in elderly, malnourished, debilitated patients, or in patients suffering from conditions which may delay wound healing.

Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. As an absorbable suture, **DAMACRYL VET** suture may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. Discard opened packages and unused sutures.

Do not re-use. Do not re-sterilize.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

DAMACRYL VET suture, which are treated to enhance handling characteristics, require the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon.

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Under some circumstances, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon.

consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur.

To avoid damaging needle points and swage areas, grasp the needle in the area one-third (1/3) to one-half (1-2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.

Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

Discard used needles in "sharps" container. Avoid prolonged exposure to elevated temperatures.

Don't use after expiry date.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occur failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occur, and transitory local irritation at the wound site.

Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. As an absorbable suture, **DAMACRYL VET** suture may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. Discard opened packages and unused sutures.

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